

GenTegraDNA™ in the Forensics Lab




GenTegraDNA for Long Term storage

All forensic cases require DNA samples to remain intact for future testing. These studies were conducted to confirm GenTegraDNA prevents degradation over long-term storage at ambient temperatures and samples remain intact for downstream applications and analysis. Freezing is the traditional gold standard for long-term storage of extracted DNA. Long term freezer storage of extracted DNA samples comes with inherent challenges, including further loss through repeated freeze-thaw cycles, adherence to the tubes, evaporation, and degradation. When DNA samples are dried in GenTegraDNA they are ready for decades of storage at room temperature or for shipping at ambient. In the dry state the GenTegraDNA coats the DNA with a thick protective layer that prevents damage to the DNA. Recovering your DNA sample after storage or shipping is easy. Simply add water to recover 100% of the DNA and the DNA is immediately ready for downstream analysis.

The Solution for Long-Term DNA Storage: Active Chemical Protection™

GenTegraDNA is a robust and reliable dry storage protection product that utilizes GenTegra's proven, patented Active Chemical Protection™ to protect DNA from hydrolysis and exposure to oxidation, and allows for higher recovery rates. GenTegraDNA provides protection even in the presence of extreme temperature conditions (-80 °C to +72 °C). Additional benefit, samples stored or shipped in GenTegraDNA do not require any specialty humidity control



GenTegraDNA does not inhibit qPCR amplification.


C_t values showed normal amplification with GenTegraDNA-protected samples compared to unprotected samples without GenTegraDNA.

Avoid hassles and risks of freezer storage

Everyone knows that freezer storage works and keeps samples stable, however, we have gotten accustomed to overlooking the daily hassles of using freezer storage. First and foremost is the risk that at any moment the freezer can fail due to mechanical issues or because the power grid fails. Also, accessing a specific sample is not straight forward as over time they all become covered in frost making identification of the desired sample difficult. You not only reduce risks and hassles, but you are also adopting a GREEN technology by using GenTegraDNA rather than freezer storage.

Downstream applications supported:

- Short Tandem Repeats, **STR**
- Next Generation Sequencing, **NGS**
- Quantitative PCR, **qPCR**
- Sample concentration




DNA samples stored on GenTegraDNA show no degradation.

Accelerated stability studies show DNA sample protection with no visible degradation after the equivalent of 20 years storage at ambient temperature.

Easy concentration of forensic samples

When using GenTegraDNA, the drying and rehydration steps make altering the sample concentration easy, without any additional steps. If the original sample of 75µL is dried

and then rehydrated in 25µL, the sample is then 3 times its original concentration. No sample loss, no guessing of how much water you have removed, and no inhibition of downstream applications!



Samples with low starting inputs were analyzed for percent recovery following GenTegra-protected and unprotected dry-down concentration methods. Seven samples (orange) were unprotected during dry-down and compared to four samples (green) protected with GenTegraDNA for the dry-down concentration process. Overall, the GenTegraDNA samples had a higher percentage of recovery, on average >200% compared to unprotected DNA samples following reconstitution.



GenTegra DNA™

Product Specification	Description
Total DNA application amount	0.00 µg – ≤20 µg
Sample application volume	20-100 µL typical range
Recovery volume	Equals application volume (20 – 100 µL of molecular biology water) Concentrating sample supported by using less than original volume
Stability for transport	Tolerance for extreme temperatures and extreme temperature shifts (-80 °C to 76 °C) Exceeds Military specifications (-60 °C to 71°C) Exceeds Federal Express® specifications (-51 °C to 60 °C)
Shelf life	3 years (prior to use)
Drying method and time	FastDryer™: Overnight SpeedVac®: 2 – 4 hours, depending on volume/type of SpeedVac Under Biosafety Hood: 14 hours
Recovery	>95%